

Date: _____ PO#: _____ Account #: _____

Doctor: _____

Address: _____

City: _____ State / Zip: _____ Phone: _____

PATIENT NAME: _____ Med Rec #: _____

Dx: _____ Dr. Email: _____

For PROLAB
office use only

P³ Functional and Diabetic Prefabs

Prices and Policies Subject to Change

P ³ Functional Adult		Price per pair	6W	7W	8W/ 6M	9W/ 7M	10W/ 8M	11W/ 9M	12W/ 10M	13W/ 11M	12M	13M
Unposted	P ³ Unposted No Topcover											
	P ³ Unposted TC Tricolor EVA Topcover											
Posted	P ³ Posted No Topcover											
	P ³ Posted TC Tricolor EVA Topcover											
	P ³ Sport High-Rebound EVA/ Microsuede Topcover											
	P ³ Low-Impact Plastazote/Poron Topcover											
Additions	Additions For P ³ Functional prefabs with topcovers only		<input type="checkbox"/> Metatarsal Pads <input type="checkbox"/> Reverse Morton's Extensions									
Diabetic	Accommodative Diabetic											13 14
Kiddythotics	Kiddythotics®		A purple	B yellow	C orange	D green	E red	F blue	G grey	Sizer Set - Includes single Kiddythotic of each size		

Matrix Prefab AFO

	Price (per AFO)	S	M	L	Max Only XL
<input type="checkbox"/> Matrix AFO					
<input type="checkbox"/> Matrix Max AFO					
<input type="checkbox"/> Right					
<input type="checkbox"/> Left					

SHIPPING AND BILLING INFORMATION

Ship via: FedEx Ground (economy) Standard Shipping

Payment: VISA MasterCard AMEX Invoice

Card # _____

Expiration ____ / ____ CVV # _____

Signature _____

PROLAB PREFABRICATED ORTHOSES

Specifications and Policies

P³ Functional Adult Prefabs

All devices are made of a semirigid polypropylene shell and functionally corrected with a 2mm medial heel skive, 14mm heel cup, and 3° valgus correction.

P³ Unposted - Functional prefab with 0/0 spot grind

P³ Unposted TC - Functional prefab with 0/0 spot grind and full-length tricolor EVA topcover

P³ Posted - Functional prefab with 4/4 polypropylene rearfoot post

P³ Posted TC - Functional prefab with 4/4 polypropylene rearfoot post and full-length tricolor EVA topcover

P³ Sport - Functional prefab with 4/4 polypropylene rearfoot post and full-length, synthetic suede topcover padded with high-rebound EVA

P³ Low-Impact - Functional prefab with 4/4 polypropylene rearfoot post and full-length, 1/8" plastazote and 1/16" poron topcover

Metatarsal pads and reverse Morton's extensions may be added to customize P³ prefabs with topcovers. (No returns on customized prefabs)

Accommodative Diabetic

Accommodative, heat-moldable prefab with P-Cell[®] top layer and Microcel Puff[®] bottom layer

P³ Functional Kiddythotics[®]

Functional prefab made for children, with: a rigid polypropylene shell, 4/4 polypropylene rearfoot post, 4mm medial heel skive, deep heel cup, and medial flange. *Sizes are color-coded.* The **Kiddythotics Sizer Set** includes a single Kiddythotic in each of the seven sizes for easy sizing.

Matrix/Matrix Max AFO (code L1932)

Prefabricated AFO for non-spastic drop foot. Highly customizable with a height-adjustable and heat-moldable anterior shell, elongated lateral strut, grindable footplate, and two washable pads. The Matrix AFO, made of dynamic carbon and glass fiber, is designed for sedentary or smaller patients. The Matrix Max, made of rigid carbon and glass fiber, is designed for active patients. Either a prefabricated or custom, topcovered foot orthosis should be worn with the AFO.

Size	Tibia Height	Foot Length	Men's Shoe Size	Women's Shoe Size
Small	14-16.5"	Up to 25cm	Up to 7	Up to 9
Medium	15-17"	24-27cm	7 to 9	9 to 10.5
Large	16-18"	25-29cm	9 to 10	10.5 to 12
X-Large (Matrix Max only)	16-19"	26-31cm	10 and Up	11.5 to 14

Policies

Prefab pricing is not subject to any discounts. Returns for unused prefabs must be made within 30 days of shipping and authorized in advance by a ProLab Client Services Representative. No returns on customized prefabs. A restocking fee is charged for each returned order. Shipping charges are the responsibility of the customer.