

Date: \_\_\_\_\_ PO#: \_\_\_\_\_ Account #: \_\_\_\_\_

Doctor: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State / Zip: \_\_\_\_\_ Phone: \_\_\_\_\_

PATIENT NAME: \_\_\_\_\_ Med Rec #: \_\_\_\_\_

Dx: \_\_\_\_\_ Dr. Email: \_\_\_\_\_

For PROLAB  
office use only

## P<sup>3</sup> Functional and Diabetic Prefabs - Special Offer

Prices and Policies Subject to Change

P <sup>3</sup> Functional Adult		Price per pair	6W	7W	8W/ 6M	9W/ 7M	10W/ 8M	11W/ 9M	12W/ 10M	13W/ 11M	12M	13M
Unposted	P <sup>3</sup> Unposted No Topcover	<del>\$14.50</del> \$13.05										
	P <sup>3</sup> Unposted TC Tricolor EVA Topcover	<del>\$20.95</del> \$18.86										
Posted	P <sup>3</sup> Posted No Topcover	<del>\$19.50</del> \$17.55										
	P <sup>3</sup> Posted TC Tricolor EVA Topcover	<del>\$25.95</del> \$23.36										
	P <sup>3</sup> Sport High-rebound EVA/ Microsuede Topcover	<del>\$29.95</del> \$26.96										
	P <sup>3</sup> Low-Impact Plastazote/Poron Topcover	<del>\$29.95</del> \$26.96										
Additions	Additions For P <sup>3</sup> Functional prefabs with topcovers only	<del>+\$10.00</del> +\$9.00 pair	<input type="checkbox"/> Metatarsal Pads <input type="checkbox"/> Reverse Morton's Extensions									
Diabetic	Accommodative Diabetic	<del>\$19.50</del> \$17.55										13   14   
Kiddythotics	Kiddythotics®	<del>\$29.50</del>	A purple	B yellow	C orange	D green	E red	F blue	G grey	Sizer Set - \$29.50 Includes single Kiddythotic of each size		
		\$26.55										

### Limited-Time Offer (June 20 - July 31, 2022)

Order prefabs between June 20<sup>th</sup> - July 31<sup>st</sup> and you are eligible to receive a complimentary pair of prefabs. Please see reverse side for details. Some restrictions apply. Make your complimentary selection below.

Style (circle 1): Unposted    Unposted TC

Posted    Posted TC

Size: \_\_\_\_\_ women's / men's

### SHIPPING AND BILLING INFORMATION

Ship via:     UPS Ground       UPS 2-Day

Payment:     VISA     MasterCard     AMEX     Invoice

Card # \_\_\_\_\_

Expiration \_\_\_\_ / \_\_\_\_      CVV # \_\_\_\_\_

Signature \_\_\_\_\_

# PROLAB PREFABRICATED ORTHOSES

## Specifications and Policies

### P<sup>3</sup> Functional Adult Prefabs

All devices are made of a semirigid polypropylene shell and functionally corrected with a 2mm medial heel skive, 14mm heel cup, and 3° valgus correction.

**P<sup>3</sup> Unposted** - Functional prefab with 0/0 spot grind

**P<sup>3</sup> Unposted TC** - Functional prefab with 0/0 spot grind and full-length tricolor EVA topcover

**P<sup>3</sup> Posted** - Functional prefab with 4/4 polypropylene rearfoot post

**P<sup>3</sup> Posted TC** - Functional prefab with 4/4 polypropylene rearfoot post and full-length tricolor EVA topcover

**P<sup>3</sup> Sport** - Functional prefab with 4/4 polypropylene rearfoot post and full-length, synthetic suede topcover padded with high-rebound EVA

**P<sup>3</sup> Low-Impact** - Functional prefab with 4/4 polypropylene rearfoot post and full-length, 1/8" plastazote and 1/16" poron topcover

*Metatarsal pads and reverse Morton's extensions may be added to customize P<sup>3</sup> prefabs with topcovers.  
(No returns on customized prefabs)*

---

### Accommodative Diabetic

Accommodative, heat-moldable prefab with P-Cell<sup>®</sup> top layer and Microcel Puff<sup>®</sup> bottom layer

---

### P<sup>3</sup> Functional Kiddythotics<sup>®</sup>

Functional prefab made for children, with: a rigid polypropylene shell, 4/4 polypropylene rearfoot post, 4mm medial heel skive, deep heel cup, and medial flange. *Sizes are color-coded.* The **Kiddythotics Sizer Set** includes a single Kiddythotic in each of the seven sizes for easy sizing.

---

### Limited-Time Offer (June 20-July 31, 2022)

When you purchase P<sup>3</sup> Prefabs between June 20-July 31, 2022, you are eligible to receive one complimentary pair, per account, of the following styles: P<sup>3</sup> Unposted, Unposted TC, Posted, or Posted TC. You must complete the complimentary pair information on the reverse side of this order form at the same time as the paid order. The paid order and complimentary pair ship together to the same address. The 10% discount applies to all prefab orders received between June 20-July 31, 2022.

---

### Policies

Prefab pricing is not subject to any additional discounts. Returns for unused prefabs must be made within 30 days of shipping and authorized in advance by a ProLab Client Services Representative. No returns on customized prefabs. A restocking fee is charged for each returned order. Shipping charges are the responsibility of the customer.